Chemie, Farben & Kosmetik

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Chemicals, Paints & Cosmetics

APPLICATION IN DETAIL

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It is never too late for better quality!


Optimisation of a unit for the manufacture of fine chemicals


Options for the upgrading of filtration in existing units taking the example of a multiple bag filter housing. It is possible to integrate the advantages of depth filter cartridges and pleated media into existing bag filter housings.

In order to ensure purity of the product, bag filter housings are used as process filters for fine chemicals for the separation of solids from liquids at all kinds of points. Due to their simple operation and their favourable cost-throughput capacity ratio, bag filter housings are therefore always the first choice in the conceptual designing of the process and the subsequent planning of the plant. Filtration must be assured and the requirements for quality assurance satisfied when the unit is started up and during subsequent operation. If the unit had been well designed, filtration with bag filters initially fulfils the desired purpose as well.

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
But often the unit no longer meets the increased requirements. In the case of a manufacturer of cosmetic products an existing process filtration section for hair shampoo was to be retrospectively optimised. The requirements for the purity of the product and hence of the unit had increased. In addition, the service life of the bag filters used up until then was unsatisfactory. The solution to the problem that suggested itself initially was to completely replace the existing filtration system with bag filters by one using filter cartridges. Such a move entails a rebuilding of the unit and is usually very costly. The addition of additional filtration steps was also discussed as a further alternative.

But it is often very difficult to implement the integration of additional filtration steps into an existing unit! Additional difficulties can arise if the conversion to multi-stage filtration with additional, new filter housings to hold filter cartridges – with better separation rates – cannot be done because there is no space or if the existing pumping capacity is no longer adequate. In such a case all the piping must be revised for the conversion.

And Wolftechnik gave some thought to the matter at precisely this point. The question being asked was: How can the advantages of depth filter cartridges and pleated media, which need to be installed within the cartridge filter housing due to their shape, also be integrated into the existing bag filter housing so that the existing unit does not need to be rebuilt?

The answer to this question led to a series of new products.

On the one hand, in the form of pleated bag filters with a larger filter area so as to lengthen the service life. On the other hand, the innovative LUPUS® II filter media from Wolftechnik found their origin here. Using both solutions, it was possible to lengthen the service life and improve the quality of the filtration of hair shampoo at the manufacturer of cosmetic products without changing anything in the unit or the existing bag filter housing!

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There was a test phase at the beginning.


Wolftechnik made the test equipment available for this purpose. Wolftechnik always has test housings available to hand. Filter media that were not needed and were still in the original packing are taken back. This means that the costs for the unit tests remain reasonable for the user. Because any form of plant conversion, also for test purposes, is a risk for the user in view of the conversion costs. It is also possible that the unit needs to be restored to its original state afterwards. Because it is only possible to definitively assess the effects on the process and the success of the action through operation in practice.

In the first step, tests were made with various kinds of filter inserts in the existing unit and the results determined regarding the service life and filtrate quality. Since no ferritic particles need to be separated out in the process filtration of hair shampoo, a test with an additional magnetic rod insert was not needed.

The most important criterion in the evaluation of the tests was an improvement in quality. In addition, the costs of the filtration in relation to the throughput capacity were assessed. LUPUS® II Multi-Layer pleated media gave the best results here (see the graphic „Service life and filtrate quality“).

The unit was successfully changed over to LUPUS® II Multi-Layer pleated media, which were integrated into the existing bag filter housing. As a result of the change-over, the manufacturer of the cosmetic products was able to meet his own more stringent requirements for quality and this without conversion costs, and the service life of the filter inserts between the individual filter change cycles was markedly longer.

The result can be transferred to all kinds of products in the chemical industry (examples: cleaning agents, body wash and shower gel, cosmetics, hygiene products, plasticisers for further processing in the plastics industry, food technology, paints and varnishes).

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Improvements through innovative products:

Due to the use of pleated bag filters of type WFB-1AP, the pleating results in an increase in service life according to the larger filtration area that is made available. At the same time, the retention rates were improved by using two-layer bag filters with internal pleating of type WFB-2IP.

A qualitative upgrade is made possible through the use of LUPUS® II filter media. Lupus®II Filter media for bag filter housings offer entirely new perspectives concerning filtrate quality and filter service life compared to conventional filter bags. The filter media can be used in the existing bag filter housings of sizes 1 and 2 without removing the filter sieve. Depending on the filtration application, different types of filter media are available for the Lupus®II filter media. They can optionally be designed as depth filter media, pleated filter media or multilayer-pleated media.

A magnetic rod can be used to separate out ferritic metal particles in addition to the pleated bag filters or the LUPUS® II filter media. Which leads here to a further lengthening of the service life.


WFB-2IP-Filterbeutel, innen doppellagig gefaltet
WFB-2IP-bag filters with double-layer internal pleating

WFB-1AP-Filterbeutel, innen einlagig gefaltet
WFB-1AP-bag filters with single-layer internal pleating

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Decision criteria for individual measures in comparison

Action	Conversion of unit	Costs of conversion	Filtration costs related to throughput	Test costs	Service life:	Filtrate quality
Additional cartridge filter housing as post-filter	Yes	High costs	Additional costs for 1 set filter cartridges	High costs, test only possibly after conversion	Remains the same	Improvement through post-filtration
Additional prefilter as strainer filter	Yes	High costs	One-time procurement costs. Additional costs for cleaning of strainers	High costs, test only possibly after conversion	Longer working life due to reducing the load on the bag filters due to the prefilter	No improvement
Replacement of the bag filter housing by a new cartridge filter housing	Yes	Moderate costs	Costs for 1 set of filter cartridges higher than for 1 set bag filters	High costs, test only possibly after conversion	Low, cartridge filters have a better retention rate than bag filters	Marked improvement
Replacement of the bag filters with pleated bag filters type WFB-1AP	No	No costs	Low additional costs	Low costs because possible without conversion	Markedly longer working life	No improvement
Replacement of the bag filters with pleated bag filters type WFB-2IP	No	No costs	Low additional costs	Low costs because possible without conversion	Markedly longer working life	Improvement due to 2-layer filtration
Replacement of the bag filters with LUPUS® II filter media	No	No costs	Costs for 1 set of LUPUS® II higher than for 1 set bag filters	Low costs because possible without conversion	Markedly longer working life	Marked improvement
Additional inclusion of a magnetic rod	No	No costs	One-time procurement costs	Low costs because possible without conversion	Markedly longer working life if there are ferritic particles	No improvement